

SK302D-32

THREE-FUNCTION LUXURIOUS ELECTRIC CARE BED



This electric medical bed has 3 functions of elevation of users' backs/legs/whole bed. The superb ergonomic design provides users greater comfort and safety and enhances efficiency for nurses and cares in operation and care services. Using toxic-free green materials and techniques, it has been widely applied in hospitals, private clinics and other medical institutions all over the world.

ADVANTAGEOUS CONFIGURATION:

- ✓ Upgraded medical motor system with 3 advanced motors, 1 control box and 1 controller with superb stability
- ✓ Permeable, smooth bed surface with holes for bedsores prevention (by advanced piercing)
- ✓ Durable, noiseless, economical polyurethane castors with impact resistance and rust-proof capabilities (castor surface with durable TPR materials)
- ✓ Low-position stability support system arranged under the bed, providing maximized support to different users (maximum 150kg)
- ✓ Rust-proof, durable, flexible aluminum guardrails by electroplating (with adjustable screw mechanism)
- ✓ Economical, deformation resistant P.P. bedside boards (ABS engineering plastic bedside boards are optional)
- ✓ Carer friendly urine hooks and stainless steel drip stand
- ✓ Can be equipped with back-up power (UPS), quick fold-down system (for CPR), nurse control system (ACP), mattress and food tables by customization
- ✓ Size: 2050x1000x460/720mm

*Weight: 100kg

For more information

+90 542 159 0407

Email: info@skymedtrade.com

