

**SK306D-53**

# FOUR-FUNCTION LUXURIOUS ELECTRIC CARE BED



**Size: 2160\*1000\*520/780 mm**

## **STANDARD CONFIGURATION:**

- ✓ Advanced medical motor system (3 motors, 1 control box, 1 handset)
- ✓ Low position stable support system.
- ✓ Central control castor brake system.
- ✓ Luxurious four folded guardrails.
- ✓ Bed board with flexible connection
- ✓ P.P head and foot board.
- ✓ Urine hooks
- ✓ S.S drip stand
- ✓ N.W: 125 KG

## **OPTIONAL CONFIGURATION:**

Mattress; food table ; back-up power system (UPS); quick fold-down system (CPR); nurse control system (ACP); ABS head and foot board; aluminum guardrails; S.S guardrails; S.S central-controlled castor brake system; back ICU; ABS bed board.

For more information

**+90 542 159 0407**

Email: [info@skymedtrade.com](mailto:info@skymedtrade.com)

