

SK402D-32

FIVE-FUNCTION LUXURIOUS ELECTRIC CARE BED



Size: 2135*1000*480/740 mm

STANDARD CONFIGURATION:

- ✓ Advanced medical motor system (4 motors, 1 control box, 1 handset)
- ✓ Low position stable support system.
- ✓ Double sided castor.
- ✓ aluminum gurdrails
- ✓ hole-punched Bed board
- ✓ P.P head and foot board.
- ✓ Urine hooks
- ✓ S.S drip stand
- ✓ N.W: 110 KG

OPTIONAL CONFIGURATION:

Mattress; food table ; back-up power system (UPS); quick fold-down system (CPR); nurse control system (ACP); ABS head and foot board.

For more information

+90 542 159 0407

Email: info@skymedtrade.com

