

SK691S

STEEL COMMODE WHEELCHAIR



Steel frame; Purplish red powder coating ; Detachable armrest; PVC seat; Detachable footrest; Black plastic footplate.

ADVANTAGES:

- ✓ Economical foldable commode wheelchair used in different places
- ✓ Stainless, aging-resistant steel frame by purplish red powder coating
- ✓ Ergonomic design based on the psysiological curvature of human waists providing users maximized comfort
- ✓ Foldable wheelchair for easy transporting and carrying in less space
- ✓ Ergonomic detachable armrests and footrests saving package size and space
- ✓ Durable and comfortable PVC cushions and backrest
- ✓ Durable, shock absorptive, maintenance-free 5"PVC wheels
- ✓ Economical durable black plastic footplates

MODEL AND SIZES

53	/	43	5"(PVC)	5"(PVC)	51	92	91	44	43	100	55*53*52	12.5	KG

For more information

+90 542 159 0407

Email: info@skymedtrade.com

