

SK809F

STANDARD STEEL WHEELCHAIR WITH FOLDABLE BACKREST AND FOOTREST
















Steel Frame; Chrome ; Fixed PVC armrest; PVC seat; Detachable footrest; Black aluminum footplate and plastic handrim

ADVANTAGES:

- ✓ Armrests and calf pads are upholstered and padded
- ✓ Quick release axles allow tool-free removal of rear wheelsDual axle hemi-height adjustable; Adjustable height back; Tough nylon upholstery

MODEL AND SIZES

												
86,5	60,5	46,41	Solid	8"(PVC)	47	89	104	39,5	42	100	75x20x79	14,7

For more information

+90 542 159 0407

Email: info@skymedtrade.com

