

# SK908AM

## STEEL WHEELCHAIR WITH FOLDABLE BACKREST



Steel Frame; Chrome; Flip up, PU armrest; Nylon seat; Flip up footrest; Black plastic footplate and plastic handrim. **(With optional Solid MAG Rear Wheel)**

### ADVANTAGES:

- ✓ High strength, heavy duty, stainless chromed steel wheelchair
- ✓ Aging-resistant steel frame with fast color
- ✓ Ergonomic design based on the physiological curvature of human waists providing users maximized comfort
- ✓ Detachable flip-up PVC armrests combined with foldable backrest saving package size and providing users maximized comfort
- ✓ Economical black plastic footplates
- ✓ Toggle-bracing brakes enabling users to stop the wheelchair rapidly, easily and safely
- ✓ Foldable wheelchair for easy transporting and carrying in less space
- ✓ Durable flame retardant nylon cushions
- ✓ Durable, shock absorptive, maintenance-free 8" PVC front wheels
- ✓ Durable, shock absorptive, maintenance-free 24" solid rear wheels

### MODEL AND SIZES

61-56	25	46(41)	24"(Solid)	8" (PVC)	49	94	105	40	45	100	80*29*76	175		

For more information

**+90 542 159 0407**

Email: [info@skymedtrade.com](mailto:info@skymedtrade.com)

