

SK909BM

STEEL WHEELCHAIR



Steel Frame; Silver Powder Coating; Flip up, PU armrest; Nylon seat; Detachable footrest; Black plastic footplate and plastic handrim. (Solid Mag Rear Wheel)

ADVANTAGES:

- ✓ High quality, comfortable, breathable, Nylon seat
- ✓ Steel Wheelchair, flexible, foldable design
- ✓ Detachable footrest; Black plastic footplate and plastic handrim design
- ✓ Flip up, PU armrest , relieve hand pressure.
- ✓ High quality rear wheel double layer design
- ✓ Reliable handbrakes, User-friendly braking system

MODEL AND SIZES

62-57-59	30.5	46-41-43	24"(PU)	8" (PVC)	47.5	96.5	107	40	45	110	80*26*95	18.6		

For more information

+90 542 159 0407

Email: info@skymedtrade.com

MOVEX
Enjoy Outside the Vision

SKYMED
Integrated Solutions