

# SK9144L

## ROLLATOR



Aluminum frame with 6" PVC castors

### ADVANTAGES:

- ✓ Purplish red powder coated rollator in European style
- ✓ Aluminum lightweight frame in foldable construction
- ✓ Hand brakes providing safety
- ✓ Seat and steel basket for rest and storage of stuff
- ✓ Auxiliary wheel in addition to four durable castors
- ✓ Castors with foot brakes
- ✓ Size: 76\*64\*81(-92)cm
- ✓ Inside Dimension: 62\*20\*84
- ✓ Max Weight Support: 75kg
- ✓ N.W.: 6kg

### MODEL AND SIZES

Length (cm)	Width (cm)	Height (cm)			PCS CN		CNTS 40 FCL
76	64	81-92	75	62*20*84	1	6	590

For more information

**+90 542 159 0407**

Email: [info@skymedtrade.com](mailto:info@skymedtrade.com)

**MOVEX**  
Enjoy Outside the Vision

**SKYMED**  
Integrated Solutions