

SK951AC

HEAVY DUTY WHEELCHAIR WITH FOLDABLE BACKREST



Steel Frame; Chrome; Flip up, PU armrest; Nylon cushion seat; Detachable & Elevating footrest; Aluminum footplate and plastic handrim

ADVANTAGES:

- ✓ Wheelchair is a large size wheelchair model specially produced for overweight users. It appeals to a wide range of users with its 150 kg carrying capacity.
- ✓ Having a foldable chassis and a semi-foldable back ensures ease of transport and storage.
- ✓ Adjustable knee-heel distance
- ✓ Footrests can be removed if desired
- ✓ It has anti-tipper wheels that prevent falling back, thus allowing the user to drive safely
- ✓ Front and rear wheels of Wheelchair are silicone infill

MODEL AND SIZES

73-68	32	56-51	24"(PU)	8" (PU)	57	95	117	48	42	125	90*35*88	21.5		

For more information

+90 542 159 0407

Email: info@skymedtrade.com

MOVEX
Enjoy Outside the Vision

SKYMED
Integrated Solutions