

SK956LAQ

ALUMINUM WHEELCHAIR



ADVANTAGES:

- ✓ With foldable backrest and quick release rear wheel
- ✓ Aluminum frame; silver liquid coating
- ✓ Flip up & height adjustable, PU armrest
- ✓ Nylon and mesh cushion(two layers) seat
- ✓ Detachable footrest
- ✓ Black plastic footplate.

MODEL AND SIZES

64-59	31	46(41)	24"(PU)	20x5(PU)	49.5	92	103	43	45	125	76*28*71	99		

For more information

+90 542 159 0407

Email: info@skymedtrade.com

MOVEX
Enjoy Outside the Vision

SKYMED
Integrated Solutions