

SK980LQ

WHEEL CHAIR STANDARD



With quick release rear wheel; aluminum frame; light blue liquid coating; flip up, PU armrest; nylon and mesh cushion(two layers);detachable footrest; black plastic footplate

ADVANTAGES:

- ✓ Aluminum children wheelchair
- ✓ Overall foldable design for easier package and transport
- ✓ Flip-up PU armrests
- ✓ Nylon and mesh cushion (two layers)
- ✓ Detachable footrests for minimum package and storage size
- ✓ pneumatic rear wheel
- ✓ Durable, shock absorptive, maintenance-free 6"PVC front wheels
- ✓ Economical black plastic footplates

MODEL AND SIZES

47	32	30	20"(PU)	6"(PVC)	48	83	96	33	38	100	75*36*86	12.8

For more information

+90 542 159 0407

Email: info@skymedtrade.com

MOVEX
Enjoy Outside the Vision

SKYMED
Integrated Solutions