

SK991

CHILDREN STANDING ASSIST SHELF



ADVANTAGES:

- ✓ Ergonomic standing aid for children with paraplegia at 3-7 years old for forward physical training
- ✓ High strength rust-proof A3 steel frame by orange powder coating:
 - Size of steel chair frame: 32*1.2
 - Size of steel movement arm: 40*10*1.0
 - Size of stainless steel adjusting bar: 30*30*1.0, 23*23*1.0
- ✓ Helpful for improving muscle atrophy and increasing muscle strength
- ✓ PU cover pillow and back cushion providing users comfort
- ✓ Chest pad and footrest cushion adjustable for different profiles of users
- ✓ Smooth, adjustable food table and footplates using high strength durable ABS materials
- ✓ Four castors with self-lock function providing safety and steadiness
- ✓ Length: 81cm
- ✓ Depth: 75cm
- ✓ Height: 84-105cm
- ✓ Front wheel size: 7.5"
- ✓ Inside dimension: 60*52*78cm
- ✓ Maximum supportive weight: 75kg

For more information

+90 542 159 0407

Email: info@skymedtrade.com

MOVEX
Enjoy Outside the Vision

SKYMED
Integrated Solutions