

SHOWER & COMMODE CHAIR

catalogue number **SK-6055D-KD**

PRODUCT CHARACTERISTICS:

- aluminum, lightweight and rust free.
- flip-up armrests , adjustable height.
- firm PE material seat and back.
- suction tipped foot to prevent slipping.
- stable seat and frame.
- removable commode pan holder and pan
- chair has backrest and armrests for support and comfort of patient's use.

USE:

used to assist patients in bathing and toileting, serving as both a shower chair and a commode chair.



Maximum capacity (kg)

100



Width (cm)

52



Height (cm)

73-85