

ALUMINUM ELEVATING WHEELCHAIR

catalogue number **SK950LAE**

PRODUCT CHARACTERISTICS:

- foldable, lightweight, aluminium frame of cross structure (double cross).
- soft nylon seat and backrest.
- one layer nylon sponge seat.
- double layer cushion .
- elevating Leg rest
- flip up armrest.
- footrest straps to secure the feet
- 24" solid rear wheels .
- solid Front wheel 6".
- PVC armrest pad.
- plastic footplate.
- fold back.
- anti-tipping wheels.

USE:

General use Adult wheelchair, to use at home and for strolls.



anti-tipping wheels



drop back handle



elevating Leg rest

MOVEX
Enjoy Outside the Vision

	16"	18"	20"
Seat width (cm)	40	45	50
Seat depth (cm)	39	39	39
Max wheelchair height (cm)	96	96	96
Max wheelchair width (cm)	59	64	70
Max wheelchair length (cm)	114	114	114
Width of folded wheelchair (cm)	31	31	31
Max capacity (kg)	140	140	150