

BABY COT

catalogue no. **SK-B**

PRODUCT CHARACTERISTICS:

Specification: L900 × W550 × H850 – 1000mm ± 20mm

The infant bassinet is made of transparent ABS, equipped with a diversion groove and ventilation holes. It is sturdy, durable, easy for observation and simple to clean.

The whole bed is integrally injection-molded from ABS, featuring bright color, high strength and long service life. The infant bassinet is tiltable within a range of 0 – 12 °, effectively preventing infants from regurgitation and choking after feeding and protecting the infant's respiratory safety.

Four deluxe silent casters are installed at the base, each with an independent brake, ensuring flexible movement, high stability, low noise and convenient mobility.

Available in green, pink and blue.

Standard Configuration: One transparent ABS infant bassinet, four deluxe silent casters.

Function Introduction:

- (1) Position adjustment: 0 – 12 ± 2 °
- (2) Height adjustment: 850 – 1000mm

