

ABS BEDSIDE CABINET

Catalogue number: **SK-BC01**

Compact design, durable structure constructed in ABS material while with light weight, all the points make this bedside locker to be an ideal option for patient ward area .

PRODUCT CHARACTERISTIC:

- whole ABS forming.
- cabinet top, drawer, tray, layer board, door and the base are made of ABS.
- towel hangers on the 2 sides of cabinet, that can be adjustable.
- dinning board is between the drawer and tabletop.
- thermos , washbasin and other commodity can be put in it.
- with four casters, easy to move.

USE:

Used beside the hospital bed for storing patients' personal items and medical supplies .



wheel diameter (cm)

5



Height (cm)

75



Width (cm)

48.5



Length (cm)

47.5




PRODUCT DETAILS

Material

The whole body is injection molded with 100% ABS plastic, which is non-toxic, strong, rust-proof, anti-corrosion, easy to



Towel hook



Towel rail is with concealed design and can be stored when not in use, same two hanging hooks are with same design for usage and storage convenience.


Sliding tray

With a hidden tray, there are grooves for placing water cups and chopsticks, which can increase the area of the table after sliding out.



ABS shoe holder

ABS shoe storage space which can put shoes in it.



Castor

The bedside table has wheels for easy movement and two brakes. It can be fixed after moving to the designated position and is durable.

