



## DRESSING TROLLEY

catalogue no. **SK-D2**

**PRODUCT CHARACTERISTICS:**

- made stainless steel 304
- two layers with two drawers
- four wheels two with brake
- with guard rails on 4 sides

**USE:**

Mainly used in hospitals and clinics for patient care procedures, especially during dressing changes. It is designed for safe and convenient transportation of medical instruments, consumables, and sterile materials. The guard rails prevent items from falling, while the drawer provides additional storage.



Height (cm)

91



Width (cm)

46



Length (cm)

60