

DRESSING TROLLEY

catalogue no. **SK-D5**

PRODUCT CHARACTERISTICS:

- made stainless steel 304
- 2 movable layers with one drawer
- four wheels two with brake
- 1 bowl, 1 cylinder waste bin, 1 bottle holder

USE:

Mainly used in hospitals and clinics for patient care procedures, especially during dressing changes. It is designed for safe and convenient transportation of medical instruments, consumables, and sterile materials.



Height (cm)

90




Width (cm)

50



Length (cm)


85



Removable Trays



waste bin



bottle holder