



# EMERGENCY TROLLEY

catalogue no. **SK-P02**

## PRODUCT CHARACTERISTICS:

- high-strength ABS boards, Top with stainless steel guard rail.
- aluminum columns with impact protection design and raised edges.
- work surface top covered with transparent soft plastic glass
- drawers: two 3" drawers: 3" H x 17.5" w x 21.5" L  
two 6" drawers : 6" H x 17.5" w x 21.5" L  
one 10" drawer: 10" H x 17.5" w x 21.5" L  
(inner dividers can be organized freely into different size units)
- 4 luxurious noiseless wheels 5" (two with brakes).
- ergonomic designed push-handle and efficient structure.
- centralized lock for all the drawers

## USE:

Used to store and transport emergency drugs and equipment for resuscitation and life-saving procedures.



Height (cm)

93



Weight (kg)

38.5



length (cm)

85



Width (cm)

52



## Power Socket

For equipment using normally length 3m

## Ceteral locking

central lock or disposable seal




## Drawer Label

Sealing side type label

## File box / utility container

For documents or other small items




## CPR board

size 54\*41.5\*2cm

## Oxygen Cylinder holder

Belts fixed on Cart



## Table top handle

Easy for pushing - Stainless steel material



## Rail

304 stainless steel rails

## IV Pole

Hidden height adjustable Iv pole  
4 hooks, dia 25mm



## 5" Braking Caster

Luxury noiseless castors  
2 with brakes

## Waste bin

Two waste bins with lid



## Extended Work surface

Hidden telescopic deputy working table  
380\*396mm

## Defibrillator Shelf

Defibrillator Shelf with multi position  
adjustments



## Sharp container

Storage for sharp needles