

# ELEVATING ADULT WHEELCHAIR

catalogue number **SK908AE**

## PRODUCT CHARACTERISTICS - BASIC:

- foldable steel frame (cross construction)
- upholstery made of black nylon
- foldable backrest
- chrome flip up armrest
- detachable footrest
- black plastic footplate
- dual axle mechanism to provide added rigidity.
- foldable wheelchair for easy transporting .
- 8"PVC front wheels
- 24"solid rear wheels
- plastic handrim

## USE:

General use wheelchair, to use at home and for strolls.



**drop back handle**



**Elevating Leg rest**

**MOVEX**  
Enjoy Outside the Vision

|                                 | 16"  | 18"  | 20"  |
|---------------------------------|------|------|------|
| Seat width (cm)                 | 40   | 45   | 51   |
| Seat depth (cm)                 | 40   | 40   | 40   |
| Max wheelchair height (cm)      | 94   | 94   | 95   |
| Max wheelchair width (cm)       | 56   | 61   | 68   |
| Max wheelchair length (cm)      | 105  | 105  | 105  |
| Width of folded wheelchair (cm) | 25   | 25   | 32   |
| wheelchair weight (kg)          | 17,5 | 18,5 | 21.2 |
| Max capacity (kg)               | 130  | 130  | 130  |