



ELEVATING CHILDREN WHEELCHAIR

catalogue number **SK980AE**

PRODUCT CHARACTERISTICS:

- foldable steel frame (cross construction)
- upholstery made of black nylon
- foldable backrest , flip up armrest
- flip up footrest
- available color Green
- black plastic footplate
- dual axle mechanism to provide added rigidity.
- foldable wheelchair for easy transporting .
- 8"PVC front wheels
- 24" spoke rear wheels
- plastic handrim.
- elevating leg rest

USE:

General use wheelchair, to use at home and for strolls.



Elevating Leg rest



drop back handle

	12"	14"
Seat width (cm)	30	35
Seat depth (cm)	33	33
Max wheelchair height (cm)	83	83
Max wheelchair width (cm)	47	52
Max wheelchair lenght (cm)	96	96
Width of folded wheelchair (cm)	32	32
Wheelchair weight (kg)	13	13
Max capacity (kg)	100	100